

TRIS AMINO™ AC TROMETHAMINE BUFFERS

Advanced Crystals
for Improved Handling
and Preparation



A Long-Standing Industry Challenge

Clumping during transport and storage is a common and challenging characteristic of conventional TRIS buffer products. There are a number of inefficiencies that can be associated with clumping of TRIS buffers, such as increased operating costs related to staff time required to de-clump product as well as longer hold times due to slower buffer solution processing. The time spent by manufacturing personnel to manually break apart aggregates of TRIS buffers and extend mixing can delay processes, reduce productivity and create operational bottlenecks.

An Innovative Solution From Advancion

To address handling challenges with TRIS buffers currently available on the market, Advancion Chemical Company developed the next-generation TRIS AMINO™ AC Advanced Crystals. TRIS AMINO AC buffers are manufactured exclusively by Advancion using proprietary process technology, and are less prone to clumping and hardness – without any change to the molecule itself. Advancion created a material solution to the clumping problem with a unique crystal morphology for TRIS AMINO AC buffers that has a larger, more rounded particle structure and does not typically form strong cakes, even when subjected to a variety of packaging, shipping and storage conditions. TRIS AMINO AC buffers may provide critical operational benefits when it comes to handling and preparation:

- Low interparticle cohesion between crystals requires less energy to break up clumps
- More spherical crystals allow for improvements in dispensing and transport through conveyors
- Lower dust emissions due to fewer fine particles
- No additives or anti-clumping agents


Why It Happens

TRIS is a hygroscopic material. It absorbs moisture, and the degree to which clumping, or “caking,” occurs is partially dependent on environmental and/or storage conditions, such as humidity, compression, temperature cycling and the duration in which these influences act on the material. Clumping is generally found as aggregates mixed in with powder but can also result in formation of solid blocks. Depending on shipping times and conditions, conventional TRIS buffers commercially available in the market can often arrive with clumping. This can also occur subsequently during storage or when drums or bags are opened and resealed during routine dispensing operations.

Larger Size, Fewer Fine Particles, Smoother Surface = Less Clumping



Conventional TRIS crystal structure



TRIS AMINO AC Advanced Crystals from Advancion

TRIS AMINO Advanced Crystals are Round and Smooth Compared to Competitive Materials Currently Available on the Market

The unique morphology of TRIS AMINO AC buffers helps reduce cake strength, resulting in less clumping and hardness than competitive materials.



Conventional TRIS crystal structure




TRIS AMINO AC Advanced Crystals from Advancion

Less Energy is Required to Break Apart Agglomerates of TRIS AMINO AC Crystals Than That Required for Competitive Materials

After 25 days exposure to high humidity and elevated temperatures, TRIS AMINO AC crystals maintain measurably lower cake strength compared to only two days' exposure for competitive products.

One-Way Analysis of Energy by Source



Product Stewardship

Through our investments in application and new product development, we have a deep understanding of where our chemistries are (and can potentially be) used to actively address global trends in Life Sciences, Beauty and Personal Care, and other essential industries to help make our lives healthier, and more sustainable and comfortable. To help ensure that Advancion products are not used in ways for which they are not intended, Advancion personnel will assist customers in dealing with environmental and product safety considerations. Before handling any of the products, obtain available product safety information including the Safety Data Sheet(s) and take the necessary steps to ensure safety of use. For assistance, product Safety Data Sheets, or other information, please visit advancionsciences.com or contact us at info@advancionsciences.com.

®/™Trademark, Advancion Corporation, registered in various countries. Notice: No freedom from infringement of any patent owned by Advancion or others is to be inferred. Because use conditions and applicable laws may differ from one location to another and may change with time, Customer is responsible for determining whether products and the information in this document are appropriate for Customer's use and for ensuring that Customer's workplace and disposal practices are in compliance with applicable laws and other government enactments. The product shown in this literature may not be available for sale and/or available in all geographies where Advancion is represented. The claims made may not have been approved for use in all countries. Advancion assumes no obligation or liability for the information in the document. References to "Advancion" or the "Company" mean the Advancion Corporation legal entity selling the products to Customer unless expressly noted. NO WARRANTIES ARE GIVEN: ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ARE EXPRESSLY EXCLUDED. September 1, 2023.